



Trouble-Free Training



ADAPTIVATION
I N C O R P O R A T E D
The World at a Touch

Pal Pads

- Can be plugged into ANY communication aid or battery-operated assistive device
- Only need 1.2 ounces of force to activate
- 3 sizes
- Easily attach objects and textures
- Great for activities
 - Adapting books
 - Attaching to cases




Flexible Switch

- Flat switch
- Adhesive-backed
- Pressure switch
- Works great in wheelchair headrests
- Flexible to the diameter of a coffee cup



SwitchBoard

- Six switches in one confined area
- Each activation area is 4" x 5"
- 2 ounces of pressure to activate
- Can be used with any device that has 1/8th-inch jacks



Taction Pads

- Not a true switch
 - Has to be used with an Adaptation product
- Activated by moisture
- Clear, adhesive-backed switch
- Can be cut!





Pneumatic Switch

- Activates with a puff of air
- Great way to make an object-based switch
- Can use squeeze toys, balls or bottles



Lex

- Three levels
- One message per level
- 20 seconds of record time for each level
- 3 volume levels
- Pairs well with TOC



Tangible Object Cards

- Encourages transitions from objects to symbols and/or text
- Great for schedules and choice making
- Child, Teen/Adult and Self-Care sets

Chipper

- Single-message device
- 20 seconds of record time
- External switch jack
- External speaker jack
- Relay jack for cause-and-effect









Sequencer

- Plays messages in sequential order
- 2 levels (independent levels!)
- 45 seconds of record time per level (short play) or 90 seconds of record time per level (long play)
- External switch jack
- Relay jack for cause-and-effect





A cow with a plow sees a sow
who takes a bow and says, "Wow!"









If You're Happy and You Know It

If you're happy and you know it,
Clap your hands.

(Clap your hands twice.)

If you're happy and you know it,
Clap your hands.

(Clap your hands twice.)

If you're happy and you know it,
And you really want to show it,
If you're happy and you know it,
Clap your hands.

(Clap your hands twice.)

If you're mad and you know it, stomp your feet. . . .
If you're excited and you know it, wave your arms. . . .
If you're sad and you know it, rub your eyes. . . .
If you're silly and you know it, shake your head. . . .





Change can make me feel scared or anxious, but I will be okay. I am learning to handle change.



My plans sometimes change. My schedule may change if people are sick, the weather is bad, or a different activity is scheduled.

Sequencing Ideas

- Tasks/Directions/Steps
- Plays
- Singing Songs
- Jokes
- Introducing people
- Ordering at a restaurant
- Requesting assistance
- Checking out library books
- Steps to make photo copies
- Listening device/memory tool

Randomizer

- Plays messages in random order
- 2 levels (independent levels!)
- 45 seconds of record time per level (short play) or 90 seconds of record time per level (long play)
- External switch jack
- Relay jack for cause-and-effect







Randomizing Ideas

- Choose reinforcer (good job, well done, great!)
- Greeting people (church, events, jobs)
- Assigning chores
- Conversation starters
- Games
- Cheering at a sporting event
- Choose your own adventure

VoicePal 8

- Up to 8 messages can be recorded and played using adaptive switches
- Delayed Activation
- 90 seconds record time (11 per message)









VoicePal 8K

- Up to 8 messages can be recorded or played using adaptive switches or the on-board keypad
- Delayed Activation
- Repeat Message Mode
- 90 seconds of record time (11 seconds each)








VoicePal Levels



- 4 Levels
- 17.5 minutes of record time
- 10 external switch jacks
- On-board keypad
- Configure for 2, 3, 5 or 10 messages
- Scanning capabilities
- Sequencing mode
- Randomizing mode







Battery Interrupter




LinkSwitch

- Used to activate battery-operated devices using switches
- Dual-channel modes
- Modes of control: Direct, Latch, Timing (seconds or minutes) and Alarm









Resources

- Recipes for Success Remixed
 - Sample recipes available at www.adaptivation.com
- Application Photo Album of Ideas
- Workshop Kit
 - (includes shipping and insurance both ways in the Continental US)
- Trial Program
 - (includes shipping and insurance one way in the Continental US)